



Product Datasheet

Phospho-STAT3 (Y705) Rabbit Monoclonal Antibody (orb548556)


Description	Phospho-STAT3 (Y705) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human STAT3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	88068 MW
Uniprot ID	P40763
Dilution Range	WB 1:5000-1:10000IHC 1:50-1:100ICC/IF 1:50-1:000IP 1:30
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Phospho-STAT3 (...)



Immunofluorescent analysis of HeLa cells...