



Product Datasheet

Phospho-GSK3 beta (Ser9) Rabbit Monoclonal Antibody (orb548617)


Description	Phospho-GSK3 beta (Ser9) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Phospho-GSK3 beta (Ser9)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	46744 MW
Uniprot ID	P49841
Dilution Range	WB 1:5000-1:10000IHC 1:50-1:200ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of GSK3 beta (phos...



IF analysis of immunocytochemical sectio...