



Product Datasheet

Phospho-ATF2 (T71) Rabbit Monoclonal Antibody (orb548618)

Description	Phospho-ATF2 (T71) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-ATF2 (T71)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	54537 MW
Uniprot ID	P15336
Dilution Range	WB 1:500-1:2000 ICC/IF 1:50-1:200 IP 1:30
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-ATF2 (T71) in SH-SY5Y cells.



ICC staining of Phospho-ATF2 (T71) in SH-SY5Y cells.