

Product Datasheet

Phospho-AMPK alpha 1 (S496) PRKAA1 Rabbit Monoclonal Antibody (orb548619)

Description Phospho-AMPK alpha 1 (S496) PRKAA1 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ICC, IF, IP, WB

Immunogen A synthesized peptide derived from human AMPK alpha 1

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 64009 MW

Uniprot ID **Q13131**

Dilution Range WB 1:1000-1:2000 ICC/IF 1:50-1:200 IP 1:50

Expiration Date 12 months from date of receipt.



Western blot analysis of Phospho-AMPK al...