


Product Datasheet

Phospho-Chk1 (S296) CHEK1 Rabbit Monoclonal Antibody (orb548625)

Description	Phospho-Chk1 (S296) CHEK1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Phospho-Chk1 (S296)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	54434 MW
Uniprot ID	014757
Dilution Range	WB 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-Chk1 (S...