

Product Datasheet

Phospho-LAT (Y191) Rabbit Monoclonal Antibody (orb548639)

Description Phospho-LAT (Y191) Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human Phospho-LAT (Y191)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 27930 MW

Uniprot ID **O43561**

Dilution Range WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:30

Expiration Date 12 months from date of receipt.



Western blot analysis of Phospho-LAT (Y191)