

Product Datasheet

Phospho-RPA2 (T21) Rabbit Monoclonal Antibody (orb548644)

Description Phospho-RPA2 (T21) Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IP, WB

Immunogen A synthesized peptide derived from human Phospho-RPA2 (T21)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 29247 MW

Uniprot ID **P15927**

Dilution Range WB 1:5000-1:10000IP 1:20

Expiration Date 12 months from date of receipt.



Western blot analysis of Phospho-RPA2 (T...