

Product Datasheet

Phospho-Tyrosine Rabbit Monoclonal Antibody (orb548666)

Description Phospho-Tyrosine Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IP, WB

Immunogen A synthesized peptide derived from Phosphotyrosine

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype Rabbit IgG

Clonality Monoclonal

Dilution Range WB 1:500-1:2000 ICC/IF 1:50-1:200 IP 1:50 FC 1:50

Expiration Date 12 months from date of receipt.



Western blot analysis of Phosphotyrosine...