

Product Datasheet

Bovine IGF2 ELISA Kit (Ready to Use) (orb548689)

Description	Bovine Insulin Like Growth Factor 2 (IGF2) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Bovine
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P07456
Sensitivity	1.44ng/mL
Entrez	281240
Expiration Date	Please enquire.