

Product Datasheet

Bovine IL15 ELISA Kit (Ready to Use) (orb548695)

Description	Bovine Interleukin 15 (IL15) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q28028
Sensitivity	6.0pg/mL
Entrez	281248
Expiration Date	Please enquire.