

Product Datasheet

Bovine TNFa ELISA Kit (Ready to Use) (orb548702)

Description	Bovine Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Bovine
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q06599
Sensitivity	32pg/mL
Entrez	280943
Expiration Date	Please enquire.