

Product Datasheet

Bovine VEGFD ELISA Kit (Ready to Use) (orb548703)

Description	Bovine Vascular Endothelial Growth Factor D (VEGFD) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	A7MBB6
Sensitivity	10.2pg/mL
Entrez	286799
Expiration Date	Please enquire.