

Product Datasheet

Canine TNFa ELISA Kit (Ready to Use) (orb548713)

Description	Canine Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Canine
Range	3.12-200pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P51742
Sensitivity	1.19pg/mL
Entrez	403922
Expiration Date	Please enquire.