

Product Datasheet

Canine VEGFA ELISA Kit (Ready to Use) (orb548715)

Description	Canine Vascular Endothelial Growth Factor A (VEGFA) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9MYV3
Sensitivity	5.4pg/mL
Entrez	403802
Expiration Date	Please enquire.