

Product Datasheet

Guinea pig BDNF ELISA Kit (Ready to Use) (orb548731)

Description	Guinea pig Brain Derived Neurotrophic Factor (BDNF) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Guinea pig
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	070183
Sensitivity	0.059ng/mL
Entrez	100730257
Expiration Date	Please enquire.