

Product Datasheet

Guinea pig CRP ELISA Kit (Ready to Use) (orb548732)

Description	Guinea pig C Reactive Protein (CRP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P49254
Sensitivity	7.4pg/mL
Entrez	100727100
Expiration Date	Please enquire.