

Product Datasheet

Guinea pig IL4 ELISA Kit (Ready to Use) (orb548746)

Description	Guinea pig Interleukin 4 (IL4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Guinea pig
Range	1.56-100pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	G0WM66
Sensitivity	0.57pg/mL
Entrez	100720403
Expiration Date	Please enquire.