

Product Datasheet

Guinea pig IL5 ELISA Kit (Ready to Use) (orb548747)

Description	Guinea pig Interleukin 5 (IL5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O08987
Sensitivity	6.4pg/mL
Entrez	100135558
Expiration Date	Please enquire.