

Product Datasheet

Guinea pig RANTES ELISA Kit (Ready to Use) (orb548750)

Description	Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) Sandwich ELISA K...
Reactivity	Guinea pig
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P97272
Sensitivity	0.059ng/mL
Entrez	100135495
Expiration Date	Please enquire.