

Product Datasheet

Guinea pig TGF beta 1 ELISA Kit (Ready to Use) (orb548751)

Description	Guinea pig Transforming Growth Factor Beta 1 (TGF beta 1) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Guinea pig
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z1Y6
Sensitivity	7.6pg/mL
Entrez	100135628
Expiration Date	Please enquire.