

Product Datasheet

Human CRISPLD2 ELISA Kit (orb548785)

Description	Human Cysteine-Rich Secretory Protein LCCL Domain-Containing 2 (CRISPLD2) Sandwich ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	1.12ng/mL
Expiration Date	Please enquire.