

Product Datasheet

Human SMURF1 ELISA Kit (orb548795)

Description	Human E3 Ubiquitin-Protein Ligase SMURF1 (SMURF1) Sandwich ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	31.3pg/mL
Expiration Date	Please enquire.