

Product Datasheet

Human EGFL7 ELISA Kit (orb548796)

Description	Human EGF Like Domain Protein, Multiple 7 (EGFL7) Sandwich ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	14.8pg/mL
Expiration Date	Please enquire.