

Product Datasheet

Human ChT ELISA Kit (orb548813)

Description	Human High Affinity Choline Transporter (ChT) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	0.057ng/mL
Expiration Date	Please enquire.