

Product Datasheet

Human IL8Ra ELISA Kit (Ready to Use) (orb548834)

Description	Human Interleukin 8 Receptor Alpha (IL8Ra) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P25024
Sensitivity	0.072ng/mL
Entrez	3577
Expiration Date	Please enquire.