

Product Datasheet

Human BSDL ELISA Kit (orb548839)

Description	Human Lipase, Bile Salt Dependent (BSDL) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	X6R868
Sensitivity	32pg/mL
Entrez	1056
Expiration Date	Please enquire.