

Product Datasheet

Human LFA3 ELISA Kit (Ready to Use) (orb548842)

Description	Human Lymphocyte Function Associated Antigen 3 (LFA3) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, cell lysates or other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	P19256
Sensitivity	5.8pg/mL
Entrez	965
Expiration Date	Please enquire.