

Product Datasheet

Human NKAP ELISA Kit (orb548858)

Description	Human NFkB Activating Protein (NKAP) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	13.8pg/mL
Expiration Date	Please enquire.