

Product Datasheet

Monkey ECF ELISA Kit (Ready to Use) (orb548917)

Description	Monkey Eosinophil Chemotactic Factor (ECF) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Monkey
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MIT7
Sensitivity	7.3pg/mL
Entrez	574218
Expiration Date	Please enquire.