

## Product Datasheet

### Monkey IFN gamma ELISA Kit (Ready to Use) (orb548926)

|                          |                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Monkey Interferon Gamma (IFN gamma) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...                                                   |
| <b>Reactivity</b>        | Monkey                                                                                                                                                    |
| <b>Range</b>             | 15.6-1000pg/mL                                                                                                                                            |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                                                              |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                                                                       |
| <b>Note</b>              | For research use only                                                                                                                                     |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.                                                   |
| <b>Uniprot ID</b>        | <b>P63309</b>                                                                                                                                             |
| <b>Sensitivity</b>       | 6.1pg/mL                                                                                                                                                  |
| <b>Entrez</b>            | <b>102128291</b>                                                                                                                                          |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                           |