

Product Datasheet

Mouse NT5E ELISA Kit (orb548945)

Description	Mouse 5'-Nucleotidase, Ecto (NT5E) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	26.1pg/mL
Expiration Date	Please enquire.