

Product Datasheet

Mouse ALDOB ELISA Kit (orb548948)

Description	Mouse Aldolase B, Fructose Bisphosphate (ALDOB) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.057ng/mL
Expiration Date	Please enquire.