

Product Datasheet

Mouse CXCL15 ELISA Kit (orb548963)

Description	Mouse Chemokine C-X-C-Motif Ligand 15 (CXCL15) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	11.8pg/mL
Expiration Date	Please enquire.