

Product Datasheet

Mouse EG-VEGF ELISA Kit (Ready to Use) (orb548975)

Description	Mouse Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) Sandwich ELISA KitReady t...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q14A28
Sensitivity	6.7pg/mL
Entrez	246691
Expiration Date	Please enquire.