

Product Datasheet

Mouse IL1 Receptor ELISA Kit (Ready to Use) (orb548999)

Description	Mouse Interleukin 1 Receptor Type I (IL1 Receptor) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13504
Sensitivity	29pg/mL
Entrez	16177
Expiration Date	Please enquire.