

Product Datasheet

Mouse IL22Ra2 ELISA Kit (Ready to Use) (orb549004)

Description	Mouse Interleukin 22 Receptor Alpha 2 (IL22Ra2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q80XF5
Sensitivity	0.119ng/mL
Entrez	237310
Expiration Date	Please enquire.