

Product Datasheet

Mouse IL8Rb ELISA Kit (Ready to Use) (orb549007)

Description	Mouse Interleukin 8 Receptor Beta (IL8Rb) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35343
Sensitivity	12.4pg/mL
Entrez	12765
Expiration Date	Please enquire.