

Product Datasheet

Mouse LDLRAP1 ELISA Kit (orb549012)

Description	Mouse Low Density Lipoprotein Receptor Adaptor Protein 1 (LDLRAP1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8C142
Sensitivity	0.067ng/mL
Entrez	100017
Expiration Date	Please enquire.