

Product Datasheet

Mouse PAWP ELISA Kit (orb549032)

Description	Mouse Postacrosomal Sheath WW Domain Binding Protein (PAWP) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	12.6pg/mL
Expiration Date	Please enquire.