

Product Datasheet

Porcine IL15 ELISA Kit (Ready to Use) (orb549070)

Description	Porcine Interleukin 15 (IL15) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Porcine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q95253
Sensitivity	12.8pg/mL
Entrez	397683
Expiration Date	Please enquire.