

Product Datasheet

Rat ALDH7A1 ELISA Kit (orb549083)

Description	Rat Aldehyde Dehydrogenase 7 Family, Member A1 (ALDH7A1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.133ng/mL
Expiration Date	Please enquire.