


Product Datasheet

Heparan Sulfate Proteoglycan 2 antibody (orb5491)

Description	Heparan Sulfate Proteoglycan 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from human Heparan Sulfate Proteoglycan 2 (3601-3700/4391 aa)
Target	HSPG2/Heparan Sulfate Proteoglycan 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	469 kDa
Uniprot ID	P98160
Expiration Date	12 months from date of receipt.



Paraffin-
embedded
mouse
embryo
tissue fi...