

Product Datasheet

Rat IL22Ra2 ELISA Kit (Ready to Use) (orb549103)

Description	Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q7TNI4
Sensitivity	0.105ng/mL
Entrez	444986
Expiration Date	Please enquire.