

Product Datasheet

Rat MYBPC3 ELISA Kit (orb549113)

Description	Rat Myosin Binding Protein C, Cardiac (MYBPC3) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	29.8pg/mL
Expiration Date	Please enquire.