

Product Datasheet

Rat PDK2 ELISA Kit (orb549121)

Description	Rat Pyruvate Dehydrogenase Kinase Isozyme 2 (PDK2) Sandwich ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	35.6pg/mL
Expiration Date	Please enquire.