

Product Datasheet

Rat SDHA ELISA Kit (orb549126)

Description	Rat Succinate Dehydrogenase Complex Subunit A (SDHA) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	31pg/mL
Expiration Date	Please enquire.