

Product Datasheet

Bovine ADP ELISA Kit (Ready to Use) (orb549148)

Description	Bovine Adiponectin (ADP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Bovine
Range	78.1-5000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q3Y5Z3
Sensitivity	32ng/mL
Entrez	282865
Expiration Date	Please enquire.