

Product Datasheet

Bovine Adrenomedullin ELISA Kit (Ready to Use) (orb549149)

Description	Bovine Adrenomedullin (Adrenomedullin) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O62827
Sensitivity	6.4pg/mL
Entrez	280713
Expiration Date	Please enquire.