

Product Datasheet

Bovine a1AGP ELISA Kit (Ready to Use) (orb549152)

Description	Bovine Alpha-1-Acid Glycoprotein (a1AGP) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Bovine
Range	15.6-1000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, urine or other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q3SZR3
Sensitivity	5.5ng/mL
Entrez	497200
Expiration Date	Please enquire.