

Product Datasheet

Bovine aLA ELISA Kit (Ready to Use) (orb549157)

Description	Bovine Alpha-Lactalbumin (aLA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Bovine
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitymilk or other biological fluids.
Sample Types	milk and other biological fluids.
Uniprot ID	P00711
Sensitivity	0.32ng/mL
Entrez	281894
Expiration Date	Please enquire.